



# Product Specification Summary

## Bacterial Culture Media

Catalog Number	40-3320-01
Product Name	LB Miller Broth (Miller's Luria Broth, High Salt); 1 Kg
Size	1 Kg
Description	LB Miller Broth (Miller's Luria Broth, High Salt); 1 kg

### LB and derivative broths

LB broth is an all purpose media for the growth of bacterial culture. LB broth contains the enzymatic digestion product of casein commonly known as peptone (also termed as Tryptone), yeast extract, and sodium chloride. LB broth is the industry standard, nutritionally rich medium for the maintenance and culture of recombinant strains of E. coli. LB can support E. coli growth OD600 2 to 3 under normal shaking incubation condition (250 rpm). Enriched media (Terrific Broth, Super Broth) can support E. coli growth to OD600 5 to 8 under normal shaking incubation conditions. Lennox formulation has the same basic composition and is low in salt content (5 g/L). It is used for selective cultivation of certain strains of E. coli that are used in salt-sensitive antibiotic selections. Miller's modification has a relatively high salt content and is suitable for selective culture of certain strains. Both of the above LB Broth formulations are also available with agar for plating. Individual formulations for the broth and agar premixed powders are given below. Several other derivative broths with high concentrations of nutrients have been formulated for higher growth and thus more plasmid and protein yields.

**General application for all enriched media** 1. E. coli growth and propagation specifically for plasmid growth followed by DNA purification.

2. Competent cell preparation for plasmid cloning and transformation.
3. Protein expression and production.
4. Phage growth and phage DNA purification.

**Specific media usage and application** 1. LB Miller is for general E. coli growth and propagation.

2. LB Lennox specifically for E. coli growth and propagation with antibiotic selection that is sensitive to high salt.
3. Terrific Broth (TB) and Super Broth (SB) for high yield protein and plasmid purification.
4. SOB (super optimal broth) for preparation of competent cells.
5. SOC (super optimal broth with catabolic repressor) for plasmid transformation and growth of competent cells.
6. 2X YT for phage propagation.

**A complete line of convenient powdered premixed media formulations are available from Gene Link and are listed below with their composition.**

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17119>



Product Manual

[http://www.genelink.com/Literature/ps/M\\_Bacterial\\_Culture\\_V2.1.pdf](http://www.genelink.com/Literature/ps/M_Bacterial_Culture_V2.1.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Agar Type A Bacterial Culture Grade; 500 g	40-3301-50	500 g
Casein Peptone(Type 1) Bacterial Culture Grade; 500 g	40-3305-50	500 g
Yeast Extract Bacterial Culture Grade; 500 g	40-3311-50	500 g
LB Miller Broth (Miller's Luria Broth, High Salt); 250 g	40-3320-25	250 g
LB Miller Broth (Miller's Luria Broth, High Salt); 500 g	40-3320-50	500 g
LB Miller Broth (Miller's Luria Broth, High Salt); 1 Kg	40-3320-01	1 Kg
LB Lennox Broth (Lennox Broth, Low Salt); 250 g	40-3321-25	250 g
LB Lennox Broth (Lennox Broth, Low Salt); 500 g	40-3321-50	500 g
LB Lennox Broth (Lennox Broth, Low Salt); 1 kg	40-3321-01	1 Kg
LB Miller Agar (Miller's Luria Broth Agar, High Sa; 250 g	40-3322-25	250 g
LB Miller Agar (Miller's Luria Broth Agar, High Sa; 500 g	40-3322-50	500 g
LB Miller Agar (Miller's Luria Broth Agar, High Salt); 1 Kg	40-3322-01	1 Kg
LB Lennox Agar (Lennox Broth Agar, Low Salt); 250 g	40-3323-25	250 g
LB Lennox Agar (Lennox Broth Agar, Low Salt); 500 g	40-3323-50	500 g
LB Lennox Agar (Lennox Broth Agar, Low Salt); 1 kg	40-3323-01	1 Kg