

Product Specification Summary

Fluorescent Molecular Probes

Catalog Number M_QPCR
Product Name QPCR_Manual
Size 1 each
Description QPCR_Manual

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17180>



Product Manual http://www.genelink.com/Literature/ps/Control_RT_Probes_PG_V3.2.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|-----------|
| SmartBase™ Fluorescent Molecular Probes | SB-FMP | 1 each |
| QPCR_Manual | M_QPCR | 1 each |
| Fluorescent Probe, 5'-Fam/3'-Quencher, 600 nmol purified yield | 26-8121-10 | 600 nmols |