



Product Specification Summary

Omni-Clean™ DNA Concentration Systems

| | |
|--|---|
| Catalog Number | 40-4026-15 |
| Product Name | G3 DNA Wash buffer 2X concentrate; 150 mL |
| Size | 150 mL |
| Description | G3 DNA Wash buffer 2X concentrate; 150 mL |
| Component/Note | G3 DNA Wash buffer 2X concentrate; 150 mL |
| Component of Omni-Clean DNA Concentration System and Omn-Mag DNA Purification System | |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17248>



Product Manual <http://www.genelink.com/Literature/ps/M40-41XX-XX V3.2.pdf>



Product MSDS http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf



Related Products

| Product | Catalog No | Size |
|--|------------|-------------------|
| Omni-Clean™ DNA Concentration System; 100 purifications | 40-4130-10 | 100 Purifications |
| Omni-Clean™ DNA Concentration System; 500 purifications | 40-4130-50 | 500 Purifications |
| Omni-Clean™ DNA Spin Column Concentration System; 100 purification | 40-4140-10 | 100 Purifications |