



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number 40-4100-05  
Product Name Omni-Mag™ DNA Purification System Kit; 500 purifications.  
Size 500 Purifications  
Description Omni-Mag DNA Purification System Kit; 500 purifications.

**Omni-Mag™ DNA Purification System Kit** An optimized kit for DNA extractions from trace samples and specifically touch DNA samples for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17368>



Product Manual [http://www.genelink.com/Literature/ps/M40-4100\\_Omni-Mag\\_V1.4.pdf](http://www.genelink.com/Literature/ps/M40-4100_Omni-Mag_V1.4.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/SDS\\_40-4100\\_Omni-Mag.pdf](http://www.genelink.com/Literature/ps/SDS_40-4100_Omni-Mag.pdf)



### Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL