



Product Specification Summary

Genetic Tools & Reagents

Catalog Number 40-4118-08

Product Name Omni-Mag™ PMP (Paramagnetic particles); 800 µL

Size 800 uL

Description Omni-Mag™ Lysis Buffer (OCC-3); 50 mL

Component/Note Omni-Mag™ PMP (Paramagnetic particles); 800 µL

Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17370>



Product Manual http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-4118_PMP_V1.2.pdf



Related Products

| Product | Catalog No | Size |
|---|------------|--------|
| Omni-Array™ mRNA Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |