



Product Specification Summary

GeneAssays™

| | |
|----------------|---|
| Catalog Number | 41-1101-01 |
| Product Name | GeneAssays Custom qPCR Probe & Primer mix, 1000 x 25uL reactions |
| Size | 1000 Reactions |
| Description | GeneAssays Custom qPCR Probe & Primer mix, 1000 x 25uL reactions |
| Component/Note | Custom probe & primer mix as per customer supplied sequence and specifications. |

Custom probe & primer mix as per customer supplied sequence and specifications.

Probe fluorophores limited to direct coupling dyes.

GeneAssays Custom qPCR Probe & Primer mix, 1000 x 25uL reactions. Supplied dried in one tube as per customer specified primer and probe concentration.

See 41-1102-01 for fluorophores that require amino to NHS conjugation.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17425>



Product Manual http://www.genelink.com/Literature/ps/GeneAssay_Custom_qPCR_Assay_PS_V3.2.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDS_nonhaz_GeneAssays.pdf



Related Products

| Product | Catalog No | Size |
|--|------------|-----------|
| Custom SPCT (Synthetic Positive Control Template); up to 130mer; 500 pmols | 26-6100-50 | 500 pmol |
| Custom SPCT (Synthetic Positive Control Template); 131-180mer; 500 pmols | 26-6101-50 | 500 pmols |
| Custom SPCT (Synthetic Positive Control Template); 131-180mer; 2 nmols | 26-6101-02 | |