



## Product Specification Summary

### Genemer™

|                |   |
|----------------|---|
| Catalog Number | 40-2036-10  |
| Product Name   | Factor VIII Genemer Pack (Hemophilia A); 10 nmols |
| Size           | 10 nmols  |
| Description    | Gene Detection System                             |
| Component/Note | Contains 27 Genemers                              |
| Component/Note | Genemer Pack                                      |

Genemer™ are optimized primer pairs for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ product contains 10 nmols of the primer pair that is sufficient for 400 regular 50 µl amplifications. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=207>



Product Manual

[http://www.genelink.com/Literature/ps/M40-2036-10\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/M40-2036-10_Ver3.1.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

| Product  | Catalog No  | Size     |
|--|-------------|----------|
| Cystic Fibrosis 7 Mutation Panel Genotyping Kit; 1 Kit | 40-2029-07K | 1 Kit    |
| Fragile X Genemer Kit; 1 Kit                           | 40-2004-11  | 1 Kit    |
| DRPLA Genemer; 10 nmols                                | 40-2042-10  | 10 nmols |
| HGH Genemer (Human Growth Hormone); 10 nmols           | 40-2024-10  | 10 nmols |
| STS Genemer (Steroid Sulfatase); 10 nmols              | 40-2023-10  | 10 nmols |
| Mouse-coxI Genemer; 10 nmols                           | 40-1502-10  | 10 nmols |
| Rat-coxI Genemer; 10 nmols                             | 40-1503-10  | 10 nmols |
| C-Hamster-Cyt.B Genemer; 10 nmols                      | 40-1504-10  | 10 nmols |
| AG Monkey-coxI Genemer; 10 nmols                       | 40-1505-10  | 10 nmols |
| Dog-coxI Genemer; 10 nmols                             | 40-1506-10  | 10 nmols |
| Cat-coxI-For1; 10 nmols                                | 40-1507-10F | 10 nmols |
| Horse-coxI Genemer; 10 nmols                           | 40-1509-10  | 10 nmols |
| Rabbit-coxI Genemer; 10 nmols                          | 40-1508-10  | 10 nmols |
| Human-coxI Genemer; 10 nmols                           | 40-1501-10  | 10 nmols |
| X alphoid repeat Genemer; 10 nmoles                    | 40-2021-10  | 10 nmols |