



Product Specification Summary

Genetic Tools & Reagents

Catalog Number	40-5025-10
Product Name	Maleic acid buffer 10X (Buffer M 10X); 100 mL
Size	100 mL
Description	Maleic acid buffer 10X (Buffer M 10X)
Component/Note	Maleic acid buffer 10X (Buffer M 10X)

Maleic acid buffer 10X (Buffer M 10X. Maleic acid buffer solution is used in non-radioactive digoxigenin based Southern blot analysis after stringent washing post hybridization with SSC/SDS based buffer solutions. Prior to anti-DIG-AP conjugate incubation the membrane is equilibrated in 1X Maleic after hybridization and then incubated in anti-DIG-AP conjugate diluted in Buffer MB composed of Maleic Acid buffer and Blocking reagent.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=4>



Product Manual http://www.genelink.com/Literature/ps/Dig_Southern_Reagents_V3.4.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit

Omni-Marker™ Universal unlabeled; 100 uL
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |
SS40-5025-10.pdf Print Date Version : December 21, 2024

40-3005-01
Page 1 of 1

