



Product Specification Summary

GeneProber™ Digoxigenin Labeled Probes

Catalog Number 40-2032-41
Product Name Kennedy Disease (SBMA) GeneProber Digoxigenin labeled Probe; 100 uL
Size 110 uL

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=655>



Product Manual <http://www.genelink.com/Literature/ps/ps40-2032-41.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|--------|
| Kennedy Disease (SBMA) GeneProber Digoxigenin labeled Probe; 100 uL | 40-2032-41 | 110 uL |
| SCA7 GeneProber Digoxigenin labeled Probe; 110 uL | 40-2041-41 | 110 uL |
| DRPLA GeneProber Digoxigenin labeled; 110 uL | 40-2042-41 | 110 uL |