



# Product Specification Summary

## Gene Detection Systems

Catalog Number	40-2004-41
Product Name	Fragile X CGG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL
Size	110 uL
Description	Fragile X CGG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL
Component/Note	CGG Triple repeat region
Component/Note	Fragile X CGG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL

Fragile X CGG repeat genotyping GeneProber™ GLFXDig1 Probe Digoxigenin labeled Fragile X CGG triple repeat spanning region digoxigenin labeled probe for non-radioactive Southern blot detection.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=75>



Product Manual <http://www.genelink.com/Literature/ps/M40-2004-41Ver6.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Factor V Leiden Genotyping Kit; 1 Kit 100 Reactions	40-2035-10K	1 Kit
Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit	40-2036-11	1 Kit

Custom Genotyping 40-2000-01  
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |  
SS40-2004-41.pdf Print Date Version : November 22, 2024 Page 1 of 1

