



# Product Specification Summary

## GeneProber™ Unlabeled Probes

Catalog Number	40-2026-39
Product Name	Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng
Size	500 ng
Description	Myotonic Dystrophy Probe
Component/Note	Unlabeled Probe-Bam HI

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=768>



Product Manual [http://www.genelink.com/Literature/ps/M40-2026-40\\_Ver5.3.pdf](http://www.genelink.com/Literature/ps/M40-2026-40_Ver5.3.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
SCA2 GeneProber unlabeled probe; 500 ng	40-2038-40	500 ng
SCA3 GeneProber unlabeled probe; 500 ng	40-2039-40	500 ng