



Product Specification Summary

Myotonic Dystrophy

Catalog Number 40-2026-15C3
Product Name Myotonic Dystrophy GScan™ Cy3 Kit
Size 1 Kit

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=822>



Product Manual http://www.genelink.com/Literature/ps/M40-2026-15_Ver6.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|--------------|--------|
| Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng | 40-2026-39 | 500 ng |
| Myotonic Dystrophy GScan™ TET Kit | 40-2026-15TT | 1 Kit |
| Myotonic Dystrophy GScan™ HEX Kit | 40-2026-15HX | 1 Kit |