



# Product Specification Summary

## Genemer™ Kits

Catalog Number	40-2027-11
Product Name	Friedreich's Ataxia Genemer™ Kit; 1 Kit
Size	1 Kit
Description	Friedreich's Ataxia
Component/Note	GAA Repeat Amplification
Component/Note	Genotyping Kit

Genemer™ are optimized primer pairs for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ product contains 10 nmols of the primer pair that is sufficient for 400 regular 50 µl amplifications. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=875>



Product Manual

<http://www.genelink.com/Literature/ps/M40-2027-11Ver5.2.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
SCA1 Genemer Kit; 1 Kit	40-2037-11	1 Kit
SCA2 Genemer Kit; 1 Kit	40-2038-11	1 Kit
SCA3 Genemer Kit; 1 Kit	40-2039-11	1 Kit
SCA6 Genemer Kit; 1 Kit	40-2040-11	1 Kit
Kennedy Disease (SBMA) Genemer Kit; 1 Kit	40-2032-11	1 Kit
SCA7 Genemer Kit; 1 Kit	40-2041-11	1 Kit
DRPLA Genemer Kit; 1 Kit	40-2042-11	1 Kit
Sickle Cell Genemer Kit; 1 Kit	40-2001-11K	1 Kit
Rh Genemer Kit; 1 Kit	40-2002-11	1 Kit
Rh EeCc (Rh Ee and Cc exon 7 specific)Genemer™ Kit; 1 Kit	40-2003-11	1 Kit
Sry Genemer Kit; 1 Kit	40-2020-11	1 Kit
X alphoid repeat Genemer™ Kit; 1 Kit	40-2021-11	1 Kit
Y alphoid repeat Genemer™ Kit; 25 uL	40-2022-11	1 Kit
STS (Steroid Sulfatase) Genemer™ Kit;	40-2023-11	1 Kit
HGH (Human Growth Hormone) Genemer™ Kit; 1 Kit	40-2024-11	1 Kit