



# Product Specification Summary

## Genemer™ Kits

Catalog Number	40-2036-11
Product Name	Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit
Size	1 Kit
Description	Factor VIII (Hemophilia) Genemer™ Kit

Genemer™ are optimized primer pairs for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ product contains 10 nmols of the primer pair that is sufficient for 400 regular 50 µl amplifications. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=889>



Product Manual [http://www.genelink.com/Literature/ps/M40-2036-10\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/M40-2036-10_Ver3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
SCA1 Genemer Kit; 1 Kit	40-2037-11	1 Kit
SCA2 Genemer Kit; 1 Kit	40-2038-11	1 Kit
SCA3 Genemer Kit; 1 Kit	40-2039-11	