



## Product Specification Summary

### Infectious Diseases

Catalog Number	60-2004-12
Product Name	Mycobacterium tuberculosis Genemer™ Mix; 200 rxns
Size	200 rxns
Description	Genemer Mix; 200 rxns.
Component/Note	Reconstitute in 200 ul
Component/Note	Contains Internal Control

The Genemer™ Primer mix product is specifically formulated to include an internal negative control primer pair and template. The amplification of an internal control negative fragment even in the absence of the positive template assures the fidelity of the amplification system. For instance, the Genemer™ Mix will amplify an internal control fragment in a 'blank' sample. Please see gel photograph for details. This Genemer™ Mix is convenient for investigators who have developed their own amplification system and require an integrated and optimized system that includes an internal negative control in each sample. The product contains specific Genemer™ pair, primer pair for negative internal control and template for negative control amplification. The quantity supplied is sufficient for 200 regular 20 ul PCR reaction. Common pathogens Genemer™ Mix are available. Please refer to our web site for a complete listing.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=970>



Product Manual

<http://www.genelink.com/Literature/ps/PS60-20XX-XX Ver 2.1.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Hepatitis B Virus (HBV) Genemer™ Amplification Kit; 1 Kit	60-2007-11KS	1 Kit
Hepatitis B Virus Genemer™; 10 nmols	60-2007-10	10 nmols
Mycobacterium tuberculosis Genemer™; 10 nmols	60-2004-10	10 nmols
Hepatitis C Virus (HCV) Genemer™ Amplification Kit; 1 Kit	60-2003-11KS	1 Kit
Hepatitis C Virus Genemer™; 10 nmols	60-2003-10	10 nmols
HIV-1 Genemer™; 10 nmols	60-2002-10	10 nmols
Hepatitis C Virus (HCV) Genemer™ Kit 100 detection; 1 Kit	60-2003-11KL	1 Kit
Hepatitis B Virus (HBV) Genemer™ Kit 100 detection; 1 Kit	60-2007-11KL	1 Kit
HIV-1 Genemer™ Detection Kit; 1 Kit	60-2002-11	1 Kit
Mycobacterium tuberculosis (MTB) Genemer™ Amplification; 1 Kit	60-2004-11KS	1 Kit
Hepatitis B Virus Genemer™ Mix; 200 rxns	60-2007-12	200 rxns
Mycobacterium tuberculosis Genemer™ Mix; 200 rxns	60-2004-12	200 rxns
Hepatitis C Virus Genemer™ Mix; 200 rxns	60-2003-12	200 rxns
HIV-1 Genemer™ Mix; 200 rxns	60-2002-12	200 rxns