



Product Specification Summary

Omni-Pure™ Viral DNA & RNA Purification Systems

| | |
|----------------|--|
| Catalog Number | 40-3720-01 |
| Product Name | Omni-Pure™ Viral DNA Purification System; 100 rxns |
| Size | 100 rxns |
| Description | Omni-Pure™ Viral DNA Purification System |

Gene Link provides a rapid purification system for extraction of viral DNA or RNA from human bodily fluids including blood. Viral DNA or RNA is captured on a special membrane and then eluted in a low volume for direct use in qualitative and quantitative amplification protocols for detection of a pathogen. The Viral RNA purification system is ideal for small volumes of human bodily fluid samples, i.e., serum, plasma and CSF. Using the easy spin column format, purification of HIV, HCV and other RNA viruses is easily accomplished in less than 30 minutes and ready for RT-PCR amplification.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=990>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3720-01.pdf>



Product MSDS http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf



Related Products

| Product | Catalog No | Size |
|--|------------|----------|
| Omni-Pure™ Viral DNA Purification System; 100 rxns | 40-3720-01 | 100 rxns |
| Omni-Pure™ Viral DNA Purification System; 500 rxns | 40-3720-05 | 500 rxns |

Omni-Pure Viral RNA Spin Column Purification System; 100 purifications
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |
SS40-3720-01.pdf Print Date Version : December 20, 2024

40-3650-01
Page 1 of 1

