



CUSTOM OLIGONUCLEOTIDE ORDER FORM FOR LONG SEQUENCES
FAX YOUR ORDER (914) 769-1193
Order by email orders@genelink.com

Date _____

SHIP TO

Name _____

Address _____

City _____ State _____ Zip _____

Phone _____ Ext. _____

Fax _____

BILL TO

Customer # _____

P.O # _____

Address _____

City _____ State _____ Zip _____

Phone _____ Ext. _____

Oligo Sequence (5' ® 3')

ID _____

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36
37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54
55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72
73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90
91	92	93	94	95	96	97	98	99	100	101	102	103	104	105	106	107	108
109	110	111	112	113	114	115	116	117	118	119	120	121	122	123	124	125	126
127	128	129	130	131	132	133	134	135	136	137	138	139	140	141	142	143	144
145	146	147	148	149	150	151	152	153	154	155	156	157	158	159	160	161	162
163	164	165	166	167	168	169	170	171	172	173	174	175	176	177	178	179	180

SCALE: 200 nmol 1 µmol **PURIFICATION:** Gel Purified Reverse Phase Cartridge*

*RPC (cartridge) purification is NOT recommended for oligos>40mer