



Product Specifications

Custom Oligo Synthesis, antisense oligos, RNA oligos, chimeric oligos, Fluorescent dyes, Affinity Ligands, Spacers & Linkers, Duplex Stabilizers, Minor bases, labeled oligos, Molecular Beacons, siRNA, phosphonates
Locked Nucleic Acids (LNA); 2'-5' linked Oligos

Oligo Modifications

For research use only. Not for use in diagnostic procedures for clinical purposes.

Other (indicate in Comments)

| | |
|--------------------------|--------|
| Category | Others |
| Modification Code | Other |
| Reference Catalog Number | 0 |
| 5 Prime | Y |
| 3 Prime | Y |
| Internal | Y |
| Molecular Weight(mw) | 349 |



select this from the drop down field when placing an order if the modification is not listed. Please describe in the comments section details about the desired modification.