



Product Specification Summary

PCR Amplification & Analysis

Catalog Number	40-3021-11
Product Name	dNTP 2mM (10X); 1.1 mL
Size	1.1 mL
Description	10 X For PCR
Component/Note	DNase, RNase Free
Component/Note	Molecular Biology Grade

Premixed dNTPs at 2mM each dATP, dCTP, dGTP, dTTP for use in all molecular biology applications. Functionally tested at Gene Link for regular PCR, Q-PCR, long range PCR, RT-PCR and cDNA synthesis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1062>



Product Manual http://www.genelink.com/Literature/ps/M40-3021-PCR_Additives_Ver5.2.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-DNA&RNA_solutions.pdf



Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	