



Product Specification Summary

Myotonic Dystrophy

Catalog Number	40-2026-15FMS
Product Name	Myotonic Dystrophy GScan™ 6-FAM Kit (20rxns); 1 Kit
Size	1 Kit
Description	Myotonic Dystrophy GScan™ 6-FAM Kit (20rxns); 1 Kit
Component/Note	6-Fam labeled primer amplification
Component/Note	DM Fluorescent Genotyping

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1076>



Product Manual http://www.genelink.com/Literature/ps/M40-2026-15_Ver6.3.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-20XX-15.pdf



Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
Myotonic Dystrophy GScan™ TET Kit	40-2026-15TT	1 Kit
Myotonic Dystrophy GScan™ HEX Kit	40-2026-15HX	1 Kit