



Product Specification Summary

RT-PCRmer™

Catalog Number 40-1004-10
Product Name beta2 microglobulin Human RT-PCRmer ; 10 nmols
Size 10 nmols

html

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=144>



Product Manual <http://www.genelink.com/Literature/ps/ps40-1004-10.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|----------|
| beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols | 40-1023-10 | 10 nmols |
| beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols | 40-1025-10 | 10 nmols |
| beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols | 40-1014-10 | 10 nmols |