

# Product Specification Summary

## BioMed Basics

# BioMed Basics

|                |  |
|----------------|--|
| Catalog Number | GL-1001-01-NC                                    |
| Product Name   | Micro tube clam shell holder 4 place; 10 Unit Pk |
| Size           | 1 Pk   |



Description Clam shell holder PVC 4 place containers for storage and secure transportation.

### Clamshell Microtube Shipping Container

PVC Clamshell Holders 4 Place & 10 Place for Storage and Shipping. Custom designed specifically for 2mL microcentrifuge tube shipping and storage. The tube slots securely clamps the tube and will not dislodge or fall off even when opened upside down. Easy clamshell lid closes firmly. Containers can also be used for convenient freezer storage.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17014>



Product Manual

[http://www.genelink.com/Literature/ps/freezer-box\\_PS.pdf](http://www.genelink.com/Literature/ps/freezer-box_PS.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

| Product  | Catalog No    | Size     |
|--|---------------|----------|
| DNA Synthesis Column Luer Lock Set; Asstd. Colors      | GL-1010-50-ST | 50 set   |
| DNA Synthesis Column Standard Luer Lock; Natural Color | GL-1011-05-NC | 50 each  |
| DNA Synthesis Column Standard Luer Lock; Blue Color    | GL-1011-05-BL | 50 each  |
| DNA Synthesis Column Standard Luer Lock; Yellow Color  | GL-1011-05-YL | 50 each  |
| DNA Synthesis Column Insert with filter; Natual Color  | GL-1012-05-NC | 50 each  |
| DNA Synthesis Column Screw Caps; Natural Color         | GL-1013-05-NC | 50 each  |
| DNA Synthesis Column Screw Caps; Green Color           | GL-1013-05-GN | 50 each  |
| DNA Synthesis Column Screw Caps; Red Color             | GL-1013-05-RD | 50 each  |
| DNA Synthesis Column Screw Caps; Blue Color            | GL-1013-05-BL | 50 each  |
| DNA Synthesis Column Screw Caps; Yellow Color          | GL-1013-05-YL | 50 each  |
| Micro tube clam shell holder 10 place; 100 Unit Pk     | GL-1002-10-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 10 Unit Pk       | GL-1001-01-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 50 Unit Pk       | GL-1001-05-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 100 Unit Pk      | GL-1001-10-NC | 1 Pk     |
| Cap Inserts for 2mL tube caps Green; 100 unit          | BB-4123-01-GN | 100 each |