


Product Specification Summary

BioMed Basics

BioMed Basics

Catalog Number BB-3100-05-CL
Product Name Freezer Storage Box 100 Place for Microcentrifuge tubes; Clear, 5 each Pk
Size 1 Pk



Description Freezer Storage Box 100 Place for 2 mL micro tubes; 5X 10 Unit Pack

Component/Note Polypropylene 100 place box

Component/Note Freezer Storage Box 100 Place with hinged lid and alphanumeric numbering for 2 mL screw cap tubes and 1.5 mL screw cap tubes. Boxes made from virgin polypropylene withstands ultra low temperature freezers

Ultralow Freezer Storage Box. 100 Place

Durable freezer boxes that withstand ultra-low freezer temperatures of -80°C and manufactured using virgin polypropylene. Redesigned with sturdy lid with 3 hinges, alphanumeric numbered for convenient cataloging. Available in a selection of color for easy color coding.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17015>



Product Manual

http://www.genelink.com/Literature/ps/freezer-box_PS.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
DNA Synthesis Column Luer Lock Set; Asstd. Colors	GL-1010-50-ST	50 set
DNA Synthesis Column Standard Luer Lock; Natural Color	GL-1011-05-NC	50 each
DNA Synthesis Column Standard Luer Lock; Blue Color	GL-1011-05-BL	50 each
DNA Synthesis Column Standard Luer Lock; Yellow Color	GL-1011-05-YL	50 each
DNA Synthesis Column Insert with filter; Natual Color	GL-1012-05-NC	50 each
DNA Synthesis Column Screw Caps; Natural Color	GL-1013-05-NC	50 each
DNA Synthesis Column Screw Caps; Green Color	GL-1013-05-GN	50 each
DNA Synthesis Column Screw Caps; Red Color	GL-1013-05-RD	50 each
DNA Synthesis Column Screw Caps; Blue Color	GL-1013-05-BL	50 each
DNA Synthesis Column Screw Caps; Yellow Color	GL-1013-05-YL	50 each
Micro tube clam shell holder 10 place; 100 Unit Pk	GL-1002-10-NC	1 Pk
Micro tube clam shell holder 4 place; 10 Unit Pk	GL-1001-01-NC	1 Pk
Micro tube clam shell holder 4 place; 50 Unit Pk	GL-1001-05-NC	1 Pk
Micro tube clam shell holder 4 place; 100 Unit Pk	GL-1001-10-NC	1 Pk
Cap Inserts for 2mL tube caps Green; 100 unit	BB-4123-01-GN	100 each