

# Product Specification Summary

## Microcentrifuge Tubes

Catalog Number BB-4110-05-GN  
Product Name Microcentrifuge tubes QuickLock™ Cap 1.5 mL Green; 500 Each  
Size 500 each



Description Microcentrifuge tubes 1.5 mL Green  
Component/Note Microcentrifuge tubes 1.5 mL Green  
Component/Note Microcentrifuge tubes 1.5 mL Green

### **1.5 mL QuickLock™ Snap Cap & 500 µL Standard Cap Microcentrifuge Tubes**

Microcentrifuge tubes are manufactured using virgin polypropylene and packaged under dust free conditions. All are RNase, DNase and pyrogen free. QuickLock™ cap design for microcentrifuge tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions. Standard microcentrifuge tubes do not have the QuickLock™ snap cap mechanism. All microcentrifuge tubes are available in a selection of color for easy color coding.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17033>



Product Manual

<http://www.genelink.com/Literature/ps/non>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

| Product   | Catalog No    | Size     |
|---|---------------|----------|
| Cap Inserts for 2mL tube caps Green; 100 unit                             | BB-4123-01-GN | 100 each |
| Cap Inserts for 2mL tube caps Red; 100 unit                               | BB-4123-01-RD | 100 each |
| Cap Inserts for 2mL tube caps Violet; 100 unit                            | BB-4123-01-VT | 100 each |
| Cap Inserts for 2mL tube caps White; 100 unit                             | BB-4123-01-WT | 100 each |
| Cap Inserts for 2mL tube caps Yellow; 100 unit                            | BB-4123-01-YL | 100 each |
| Cap Inserts for 2mL tube caps assorted colors; 100 unit                   | BB-4123-01-AS | 100 each |
| Cap Inserts for 2mL tube caps Blue; 100 unit                              | BB-4123-01-BL | 100 each |
| Cap Inserts for 2mL tube caps Orange; 100 unit                            | BB-4123-01-OR | 100 each |
| Microcentrifuge conical screw cap tube 2mL Amber without caps; 500 Each   | BB-4120-05-AM | 500 each |
| Microcentrifuge conical screw cap tube 2mL Natural without caps; 500 Each | BB-4120-05-NC | 500 each |
| Microcentrifuge conical tube 2mL with Amber screw cap; 500 Each           | BB-4121-05-AM | 500 each |
| Microcentrifuge conical tube 2mL with Blue screw cap; 500 Each            | BB-4121-05-BL | 500 each |
| Microcentrifuge conical tube 2mL with Gray screw cap; 500 Each            | BB-4121-05-GY | 500 each |
| Microcentrifuge conical tube 2mL with Green screw cap; 500 Each           | BB-4121-05-GN | 500 each |
| Microcentrifuge conical tube 2mL with Naturalscrew cap; 500 Each          | BB-4121-05-NC | 500 each |