



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-3020-05
Product Name	dNTP 10mM; 0.5 mL
Size	0.5 mL
Description	Molecular Biology Grade
Component/Note	DNase, RNase Free

Premixed dNTPs at 10mM each dATP, dCTP, dGTP, dTTP for use in all molecular biology applications. Functionally tested at Gene Link for regular PCR, Quantitative PCR, long range PCR, RT-PCR and cDNA synthesis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17089>



Product Manual [http://www.genelink.com/Literature/ps/M40-3021-PCR\\_Additives\\_Ver5.2.pdf](http://www.genelink.com/Literature/ps/M40-3021-PCR_Additives_Ver5.2.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/SDS\\_40-DNA&RNA\\_solutions.pdf](http://www.genelink.com/Literature/ps/SDS_40-DNA&RNA_solutions.pdf)



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL