



Product Specification Summary

siRNA:RNA Interference

| | |
|----------------|---|
| Catalog Number | 27-6437-01 |
| Product Name | Custom siRNA duplex Reverse Phase Purified; 100 nmols |
| Size | 100 nmols |
| Description | Custom siRNA duplex Reverse Phase Purified, 100 nmols |

Custom siRNA are synthesized at different scales of synthesis, the yield varies between these synthesis scales, individual synthesis, modifications and purification methods. Your custom siRNA duplex is accompanied by a Certificate of Analysis stating the exact yield in nmols and μg . The exact specifications are also printed on the tube label. Please refer to the Certificate of Analysis or the tube label for the exact specifications. Place your custom siRNA order online at the following link <http://www.genelink.com/newsite/allonlineorders.asp> or you can email your sequences at oligos@genelink.com

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17111>



Product Manual http://www.genelink.com/Literature/ps/siRNA_Custom_Duplex_PS_V1.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|-------|
| Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit | 27-6402-01 | 1 Kit |
| 3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each | 27-6402-02 | 1 Kit |

Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |
SS27-6437-01.pdf Print Date Version : July 17, 2024

27-6402-06
Page 1 of 1

