



## Product Specification Summary

### Sequencing Primers

Catalog Number	26-3000-56
Product Name	T7 terminator
Size	25 ug
Description	T7 terminator
Component/Note	Sequence: 5'-GCTAGTTATTGCTCAGCGGT-3'

The T7, T3 and SP6 RNA polymerase primers are complementary to the conserved class III sequence. Sequencing data will be obtained for all templates containing the Class III RNA polymerase promoter sequences. The product is supplied as a lyophilized powder, after reconstitution store at -20oC. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17175>



Product Manual [http://www.genelink.com/Literature/ps/T7\\_T3\\_seq\\_primer\\_PSV3.1.pdf](http://www.genelink.com/Literature/ps/T7_T3_seq_primer_PSV3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



#### Related Products

Product	Catalog No	Size
M13 Rev Seq; 25 ug	26-3000-52	25 ug
M13 For Seq; 25 ug	26-3000-51	25 ug
T7 Seq Long	26-3000-55	