



## Product Specification Summary

### Omni-Clean™ DNA Concentration Systems

Catalog Number	40-4026-15
Product Name	G3 DNA Wash buffer 2X concentrate; 150 mL
Size	150 mL
Description	G3 DNA Wash buffer 2X concentrate; 150 mL
Component/Note	G3 DNA Wash buffer 2X concentrate; 150 mL
Component of Omni-Clean DNA Concentration System and Omn-Mag DNA Purification System	

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17248>



Product Manual <http://www.genelink.com/Literature/ps/M40-41XX-XX V3.2.pdf>



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_40-Omni-Pure.pdf](http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf)



#### Related Products

Product	Catalog No	Size
Omni-Clean™ DNA Concentration System; 100 purifications	40-4130-10	100 Purifications
Omni-Clean™ DNA Concentration System; 500 purifications	40-4130-50	500 Purifications
Omni-Clean™ DNA Spin Column Concentration System; 100 purification	40-4140-10	100 Purifications