

Product Specification Summary

Fluorescent Molecular Probes

Catalog Number	FPRS
Product Name	Fluorescent Molecular Probes Reconstitution & Stability
Size	1 1
Description	Fluorescent Molecular Probes Reconstitution & Stability
Component/Note	Fluorescent Molecular Probes Reconstitution & Stability
Component/Note	Fluorescent Molecular Probes Reconstitution & Stability

Reconstitution & Storage

Gene Link oligos are supplied lyophilized. These are stable at room temperature for an extended period of time. The oligonucleotide should preferably be frozen upon receipt. TE buffer (10mM Tris, 1mM EDTA, pH 7.5) is recommended for dissolving the oligonucleotides; EDTA inhibits the activity of the nucleases.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17393>

Product Manual  http://www.genelink.com/Literature/ps/Fluorescent_oligo_recons_V3.2.pdf

Product MSDS  <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
SmartBase™ Fluorescent Molecular Probes	SB-FMP	1 each
QPCR_Manual	M_QPCR	1 each
Fluorescent Probe, 5'-Fam/3'-Quencher, 600 nmol purified yield	26-8121-10	600 nmols
Fluorescent Probe, 5'-Fam/3'-Quencher, 900 nmol purified yield	26-8121-15	900 nmols
Fluorescent Probe 5'-CAL Orange 560/3'-Quencher 600 nmol purified yield	26-8142-10	600 nmols
Fluorescent Probe 5'-CAL Orange 560/3'-Quencher 900 nmol purified yield	26-8142-15	900 nmols
Fluorescent Probe 5'-CAL Red 610/3'-Quencher 900 nmol purified yield	26-8141-15	900 nmols
Fluorescent Probe 5'-CAL Red 610/3'-Quencher; 600 nmol purified yield	26-8141-10	600 nmols
Fluorescent Probe, 5'-Quaser 670/Quencher. 500 nmol purified yield	26-8134-10	500 nmols
Fluorescent Probe, 5'-Quaser 670/Quencher. 750 nmol purified yield	26-8134-15	750 nmols
Fluorescent Probe, 5'-Fam/3'-Quencher; 120 nmol purified yield	26-8121-03	120 nmols
Fluorescent Probe, 5'-Quaser 670/Quencher. 80 nmol purified yield	26-8134-03	80 nmols
Fluorescent Probe, 5'-CAL Red 610/Quencher; 120 nmols purified yield	26-8141-03	120 nmols
Fluorescent Probe, 5'-CAL Orange 560/Quencher; 120 nmol yield	26-8142-03	120 nmols
Fluorescent Probe, 5'-Hex/3'-Quencher; 50 nmols purified yield	26-8122-01	50 nmols