

Product Specification Summary

Fluorescent Molecular Probes

Catalog Number 26-8124-01
Product Name Fluorescent Probe, 5'-Cy 5/3'-Quencher; 50 nmols purified yield
Size 60 nmols
Description Fluorophore-Quencher Probe, 5'-Cy 5/3'-Quencher; 50 nmols purified yield
Component/Note Fluorescent Probe, 5'-Cy 5/3'-Quencher; 50 nmols purified yield
Scan the QR Code or visit the following links
Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17415>

Product Manual http://www.genelink.com/Literature/ps/Fluorescent_oligo_recons_V3.2.pdf

Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>

Related Products

Product	Catalog No	Size
SmartBase™ Fluorescent Molecular Probes	SB-FMP	1 each
QPCR_Manual	M_QPCR	1 each
Fluorescent Probe, 5'-Fam/3'-Quencher, 600 nmol purified yield	26-8121-10	600 nmols