



# Product Specification Summary

## Sry, X & Y

Catalog Number 44-6130-01  
Product Name X alphoid TaqMan 5'FAM; 25 ug  
Size 25 ug

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=339>



Product Manual <http://www.genelink.com/Literature/ps/PS44-6430-01.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product                                | Catalog No | Size  |
|--|------------|-------|
| X alphoid TaqMan 5'FAM; 25 ug          | 44-6130-01 | 25 ug |
| X alphoid TaqMan 5' Fluorescein; 25 ug | 44-6140-01 | 25 ug |
| SRY TaqMan 5' Fluorescein; 25 ug       | 44-6440-01 | 25 ug |