



# Product Specification Summary

## GeneProber™ Unlabeled Probes

|                |   |
|----------------|---|
| Catalog Number | 40-2026-39  |
| Product Name   | Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng |
| Size           | 500 ng  |
| Description    | Myotonic Dystrophy Probe  |
| Component/Note | Unlabeled Probe-Bam HI  |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=768>



Product Manual [http://www.genelink.com/Literature/ps/M40-2026-40\\_Ver5.3.pdf](http://www.genelink.com/Literature/ps/M40-2026-40_Ver5.3.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product   | Catalog No | Size   |
|---|------------|--------|
| Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng | 40-2026-39 | 500 ng |
| SCA2 GeneProber unlabeled probe; 500 ng                             | 40-2038-40 | 500 ng |
| SCA3 GeneProber unlabeled probe; 500 ng                             | 40-2039-40 | 500 ng |