

Product Specification Summary

SCA2

Catalog Number 40-2038-11
Product Name SCA2 Genemer Kit; 1 Kit
Size 1 Kit

Genemer™ genotyping kits are optimized components for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Infomation <http://www.genelink.com/geneprodsite/product.asp?p=801>



Product Manual http://www.genelink.com/Literature/ps/M40-2038-11_V4.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
SCA2 22 repeat (19 CAG + 3 CAA) GScan & Genemer Control DNA	40-2038-01	500 ng
SCA2 GeneProber unlabeled probe; 500 ng	40-2038-40	500 ng
SCA2 Genemer Kit; 1 Kit	40-2038-11	1 Kit