



Product Specification Summary

GeneProber™ Unlabeled Probes

Catalog Number 40-2026-38
Product Name Myotonic Dystrophy GeneProber GLDM3 unlabeled Probe; 500 ng
Size 500 ng

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=870>



Product Manual http://www.genelink.com/Literature/ps/M40-2026-40_Ver5.2.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
SCA2 GeneProber unlabeled probe; 500 ng	40-2038-40	500 ng
SCA3 GeneProber unlabeled probe; 500 ng	40-2039-40	500 ng