



## Product Specification Summary

### Sequencing Primers

Catalog Number 26-3000-11  
Product Name Lambda gt10 reverse primer 24mer; 25 ug  
Size 25 ug

This lambda gt10 reverse primer is upstream of the unique EcoRI cloning site present in the phage imm434 gene region. This primer can be used to sequence inserts, and can be used with the Lambda gt10 forward primer (catalog # 26-3000-08) to amplify the insert using polymerase chain reaction. The product is supplied as a lyophilized powder, after reconstitution store at -20 C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=89>



Product Manual [http://www.genelink.com/Literature/ps/Lambda\\_seq\\_primer\\_PSV3.1.pdf](http://www.genelink.com/Literature/ps/Lambda_seq_primer_PSV3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



#### Related Products

Product	Catalog No	Size
M13 Rev Seq; 25 ug	26-3000-52	25 ug
M13 For Seq; 25 ug	26-3000-51	25 ug
T7 Seq Long	26-3000-55	25 ug