



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number 40-3071-01  
Product Name Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb  
Size 100 ul

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=910>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3062-XX.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL