

# Product Specification Summary

## Fluorescent Molecular Probes

Catalog Number	FMPPM
Product Name	Fluorescent Molecular Probes - Product Manual
Size	1 Doc
Description	Fluorescent Molecular Probes - Product Manual
Component/Note	Fluorescent Molecular Probes - Product Manual
Component/Note	Fluorescent Molecular Probes - Product Manual
Fluorescent Molecular Probes - Product Manual	

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=948>



Product Manual [http://www.genelink.com/Literature/ps/Fluorescent\\_oligo\\_recons\\_V3.2.pdf](http://www.genelink.com/Literature/ps/Fluorescent_oligo_recons_V3.2.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
SmartBase™ Fluorescent Molecular Probes	SB-FMP	1 each
QPCR_Manual	M_QPCR	1 each
Fluorescent Probe, 5'-Fam/3'-Quencher, 600 nmol purified yield	26-8121-10	600 nmols