



## Product Specification Summary

### siRNA:RNA Interference

Catalog Number	27-6415-02FL
Product Name	Rat Chromogranin-A custom siRNA duplex FL labeled; 2 nmols
Size	2 nmols
Description	Pos. Control [NM_021655]
Component/Note	Fluorescein Labeled
Component/Note	Chromogranin siRNA Duplex

A negative control should also be included when performing siRNA experiments to eliminate the possibility of nonspecific silencing effects. For this purpose, Gene Link provides a siRNA sequence that has been shown to have no effect on gene silencing. The negative control siRNA is also available labeled with fluorescein, which allows transfection efficiency to be tracked by fluorescent microscopy.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=962>



Product Manual

[http://www.genelink.com/Literature/ps/RNAi\\_Controls\\_COA\\_PS\\_V6.1.pdf](http://www.genelink.com/Literature/ps/RNAi_Controls_COA_PS_V6.1.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit
Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit	27-6402-06	1 Kit
siRNA non-silencing FL labeled; 20 nmols	27-6411-20FL	20 nmols
Human Vimentin custom siRNA duplex; 20 nmols	27-6412-20	20 nmols
Mouse Beta Tubulin control siRNA duplex; 20 nmols	27-6414-20	20 nmols
Rat Chromogranin-A control siRNA duplex; 20 nmols	27-6415-20	20 nmols
Human Vimentin control siRNA duplex FL label; 10 nmols	27-6412-10FL	10 nmols
Human Beta Actin control siRNA duplex FL; 10 nmols	27-6413-10FL	10 nmols
Rat Chromogranin-A control siRNA duplex; 10 nmols	27-6415-10FL	10 nmols
Custom siRNA duplex standard purity; 2 nmols	27-6436-20	2 nmols
Custom siRNA duplex standard purity; 10 nmols	27-6436-10	10 nmols
Custom siRNA duplex standard purity; 50 nmols	27-6436-05	50 nmols
Custom siRNA duplex standard purity; 100 nmols	27-6436-01	100 nmols
Custom siRNA duplex Reverse Phase Purified; 2 nmols	27-6437-20	2 nmols