



Product Specification Summary

Infectious Diseases

Catalog Number	60-2002-10
Product Name	HIV-1 Genemer TM , 10 nmols
Size	10 nmols
Description	HIV-1 Primer P1/P2
Component/Note	10 nmol primer pair

Optimized primer pair for PCR amplification of HBV specific DNA fragment. Specialized amplification conditions may be required for certain amplifications. Amplification of target fragment DNA requires optimization using varying amounts of the template based on its abundance. Generally less than 10 ng of template is sufficient for a successful amplification. It is a good strategy initially to amplify multiple amounts for optimization. Use known positive control DNA as one sample. The product contains 10 nmols lyophilized GenemerTM pair. The quantity supplied is sufficient for 800 regular 25 ul PCR reaction. Common pathogen GenemerTM are available. Please refer to our web site for a complete listing.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=978>



Product Manual

<http://www.genelink.com/Literature/ps/PS60-20XX-XX Ver 2.1.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Hepatitis B Virus (HBV) Genemer™ Amplification Kit; 1 Kit	60-2007-11KS	1 Kit
Hepatitis B Virus Genemer™; 10 nmols	60-2007-10	10 nmols
Mycobacterium tuberculosis Genemer™; 10 nmols	60-2004-10	10 nmols
Hepatitis C Virus (HCV) Genemer™ Amplification Kit; 1 Kit	60-2003-11KS	1 Kit
Hepatitis C Virus Genemer™; 10 nmols	60-2003-10	10 nmols
HIV-1 Genemer™; 10 nmols	60-2002-10	10 nmols
Hepatitis C Virus (HCV) Genemer™ Kit 100 detection; 1 Kit	60-2003-11KL	1 Kit
Hepatitis B Virus (HBV) Genemer™ Kit 100 detection; 1 Kit	60-2007-11KL	1 Kit
HIV-1 Genemer™ Detection Kit; 1 Kit	60-2002-11	1 Kit
Mycobacterium tuberculosis (MTB) Genemer™ Amplification; 1 Kit	60-2004-11KS	1 Kit
Hepatitis B Virus Genemer™ Mix; 200 rxns	60-2007-12	200 rxns
Mycobacterium tuberculosis Genemer™ Mix; 200 rxns	60-2004-12	200 rxns
Hepatitis C Virus Genemer™ Mix; 200 rxns	60-2003-12	200 rxns
HIV-1 Genemer™ Mix; 200 rxns	60-2002-12	200 rxns