



## Product Specification Summary

### Omni-Pure™ Viral DNA & RNA Purification Systems

Catalog Number	40-3720-01
Product Name	Omni-Pure™ Viral DNA Purification System; 100 rxns
Size	100 rxns
Description	Omni-Pure™ Viral DNA Purification System

Gene Link provides a rapid purification system for extraction of viral DNA or RNA from human bodily fluids including blood. Viral DNA or RNA is captured on a special membrane and then eluted in a low volume for direct use in qualitative and quantitative amplification protocols for detection of a pathogen. The Viral RNA purification system is ideal for small volumes of human bodily fluid samples, i.e., serum, plasma and CSF. Using the easy spin column format, purification of HIV, HCV and other RNA viruses is easily accomplished in less than 30 minutes and ready for RT-PCR amplification.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=990>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3720-01.pdf>



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_40-Omni-Pure.pdf](http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf)



#### Related Products

Product	Catalog No	Size
Omni-Pure™ Viral DNA Purification System; 100 rxns	40-3720-01	100 rxns
Omni-Pure™ Viral DNA Purification System; 500 rxns	40-3720-05	500 rxns

Omni-Pure Viral RNA Spin Column Purification System; 100 purifications  
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS40-3720-01.pdf Print Date Version : July 17, 2024

40-3650-01  
Page 1 of 1

